**ROTIXA 500 RS**

Ⓓ	<b>Bedienungsanleitung</b> .....	<b>9</b>
Ⓔ	<b>Operating Manual</b> .....	<b>37</b>
Ⓕ	<b>Mode d'emploi</b> .....	<b>64</b>
Ⓘ	<b>Istruzioni per l'uso</b> .....	<b>92</b>

Andreas Hettich GmbH & Co. KG  
Föhrenstraße 12, D-78532 Tuttlingen / Germany  
Phone +49 (0)7461 / 705-0  
Fax +49 (0)7461 / 705-1125  
info@hettichlab.com, service@hettichlab.com  
www.hettichlab.com

© 2014 by Andreas Hettich GmbH & Co. KG

All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.

Änderungen vorbehalten! , Modifications reserved! , Sous réserve de modifications ! , Con riserva di modifiche!

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We hereby declare under our sole responsibility that the designated device and its accessories, which are listed in the technical documentation for this device and whose conformity has been assessed together with the device, conform to the Medical Device Directive – Council Directive 93/42/EEC.

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The conformity evaluation process was performed in accordance with appendix II of Directive 93/42/EEC involving the following notified bodies:

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**The following additional European directives and ordinances have been applied:**

- Machinery Directive 2006/42/EU
- EMC directive 2014/30/EU
- Low voltage directive 2014/35/EU
- RoHS II Directive 2011/65/EU (without involvement of a notified body)
- Ordinance (EC) No. 1907/2006 (REACH) (without involvement of a notified body)

Applied standards:

According to the list of applied standards, which is part of the product file.

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Last validity of this declaration:

In accordance with the last date of validity of the EC certificate of conformity as issued by the above-mentioned notified body in accordance with appendix II of Directive 93/42/EEC which was valid at the time of the issuance of this declaration: 2020-06-09

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- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Kanada:** CMDR, SOR/98-282 "Medical Devices Regulations"

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## Standards and regulations which apply to this device

The device is a high-end technical product. It is subject to extensive testing and certification procedures according to the following standards and regulations in their respectively valid version:

### Electrical and mechanical safety for design and final testing:

Standard series: IEC 61010 (conform to standards of DIN EN 61010)

- IEC 61010-1 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements" (Pollution Degree 2, Excess-voltage category II)
- IEC 61010-2-010 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials" (applied to heated centrifuges only)
- IEC 61010-2-020 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-020: Particular requirements for laboratory centrifuges"
- IEC 61010-2-101 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment"

### Electromagnetic Compatibility:

- EN 61326-1 "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements"

### Risk management:

- DIN EN ISO 14971 "Application of risk management to medical devices"

### Restriction of Hazardous Substances (RoHS II):

- EN 50581 "Technical documentation for assessing electric and electronic devices with regard to the restriction of hazardous substances"

### European directives applied for conformity assessment procedures:

- Medical device directive 93/42/EEC (class 2a medical products)  
EC conformity assessment procedure according to Annex II  
"EC DECLARATION OF CONFORMITY (complete quality assurance system)".  
The certificate, Annex II, is issued by the following notified body:

MEDCERT GmbH – Notified Body CE 0482  
tel: +49 (0)40 - 22633250  
fax: +49 (0)40 - 226332515  
e-mail: info@medcert.de  
website: www.med-cert.com  
mail: Pilatuspool 2, D-20355 Hamburg; Germany

- Directive 2011/65/EU for the restriction of use of certain hazardous substances in electric and electronic devices. Carrying out the EC conformity assessment process is the sole responsibility of the manufacturer, without the involvement of a notified body.

### Applied medical device regulations outside Europe:

- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada:** CMDR, SOR/98-282 "Medical Devices Regulations"

### Certified quality management system according to

- ISO 9001 "Quality management systems – Requirements"
- ISO13485 "Medical devices - Quality management systems - Requirements for regulatory purposes"

### Environmental management system according to

- ISO 14001 "Environmental management systems - Requirements with guidance for use"

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### Restrictions relatives aux substances dangereuses (RoHS II):

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- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs : TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada:** CMDR, SOR/98-282 "Medical Devices Regulations"

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## Contents

1	How to use these operating instructions .....	39
2	Symbol meanings.....	39
3	Intended use.....	40
4	Remaining risks.....	40
5	Technical specifications .....	41
6	Notes on safety .....	43
7	Unpacking and setting up the centrifuge.....	45
8	Delivery checklist .....	45
9	Transport and storage.....	46
9.1	Transport .....	46
9.2	Storage .....	46
10	Initial operation .....	46
11	Interface (only for centrifuge with interface).....	47
12	Opening and closing the lid.....	47
12.1	Opening the lid .....	47
12.2	Closing the lid .....	47
13	Installation and removal of the rotor.....	47
14	Loading the rotor .....	48
15	Closing biosafety systems.....	49
16	Control and display devices .....	50
16.1	Status symbols .....	50
16.2	Control knob .....	50
16.3	Keys and set-up options .....	50
17	Entering centrifugation parameters .....	51
17.1	Running time.....	51
17.1.1	Time selection.....	51
17.1.2	Continuous running.....	51
17.2	Speed (RPM).....	51
17.3	Relative centrifugal force (RCF/RZB).....	51
17.4	Run-up and run-down parameters .....	51
17.4.1	Run-up step .....	51
17.4.2	Run-up time .....	51
17.4.3	Braking step .....	52
17.4.4	Run-down time.....	52
17.4.5	Brake de-energisation speed.....	52
17.5	Radius/Temperature.....	52
17.5.1	Radius.....	52
17.5.2	Temperature .....	52
17.6	Automatic buffer .....	52
18	Programming.....	52
18.1	Entering / modifying programs .....	52
18.2	Calling up programs .....	52
19	Centrifugation .....	53
19.1	Centrifugation with time selection.....	53
19.2	Continuous running .....	53

20	Changing settings during the centrifugation run.....	53
21	Integral RCF (jRCF) .....	53
22	Displaying the maximum speed of the rotor .....	53
23	Displaying the maximum RCF of the rotor .....	54
24	Emergency stop.....	54
25	Audible signal .....	54
26	Interrogating the operating hours .....	54
27	Setting the date and time.....	54
28	Centrifugation data displayed after switching on.....	54
29	Immediate display of the centrifugation data after switching on .....	55
30	Key-operated switch.....	55
31	Program linking (for centrifuges with program linking only) .....	55
31.1	Linking programs .....	55
31.2	Changing program links.....	55
31.3	Centrifugation run with program linking .....	56
31.4	Deleting program linking .....	56
32	Cooling.....	56
32.1	Standby-cooling .....	56
32.2	Pre-cooling the rotor .....	56
33	Heating .....	56
34	Relative centrifugal force (RCF).....	57
35	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm <sup>3</sup> .....	57
36	Rotor recognition .....	57
37	Emergency release.....	58
38	Maintenance and servicing.....	58
38.1	Centrifuge (housing, lid and centrifuging chamber).....	58
38.1.1	Surface cleaning and care .....	58
38.1.2	Surface disinfection.....	59
38.1.3	Removal of radioactive contaminants.....	59
38.2	Rotors and Attachments .....	59
38.2.1	Cleaning and care.....	59
38.2.2	Disinfection.....	60
38.2.3	Removal of radioactive contaminants.....	60
38.2.4	Trunnions .....	60
38.2.5	Rotors and accessories with limited service lives.....	60
38.2.6	Calculation of the number of run cycles carried out .....	60
38.3	Autoclaving .....	61
38.4	Centrifuge containers.....	61
39	Faults .....	62
40	Returning Devices .....	63
41	Disposal .....	63
42	Anhang / Appendix .....	120
42.1	Rotoren und Zubehör / Rotors and accessories.....	120

## 1 How to use these operating instructions

- Before using the centrifuge, read the operating instructions and observe them.
- The operating instructions are a part of the device. They must always be kept readily available.
- If the device is set up at a different location, the operating instructions must be provided with it.

## 2 Symbol meanings

Symbol on the device:

Attention, general hazard area.

Before using the device, make sure you read the operating instructions and observe the safety information!

Symbol in this document:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the device and in this document:

Beware of biohazard.

Symbol on the device and in this document:

Beware of hot surface.

Nonobservance of this warning can lead to material damage and personal injury.

Symbol on the device:

Equipotential: Connector (PE connector) for potential equalization (only for centrifuge with PE connector).

Symbol on the device:

Switch positions of the key-operated switch.

Symbol on the device:

Optical interface (only for device with optical interface).

Symbol in this document:

This symbol refers to important circumstances.

Symbol on the device and in this document:

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

### 3 Intended use

This device is a laboratory centrifuge which is suitable for medical and in-vitro diagnostic applications.

The therapeutic purpose with the corresponding centrifuge vessels is the separation of blood by means of centrifugation for the purpose of transfusions or autotransfusion of the won blood components.

The centrifuge is also suitable for preparing human bodily fluids by means of centrifugation for in-vitro diagnostic purposes.

Furthermore, the centrifuge is generally also suitable for separating substances or substance mixtures of non-human origin having a density of max. 1.2 kg/dm<sup>3</sup>.

The centrifuge may only be used by skilled personnel in closed laboratories.

The centrifuge is only meant for the above-named purposes.

Another use or one which goes beyond this, is considered to be non-intended. The company Andreas Hettich GmbH & Co. KG is not liable for damage resulting from this.

Observing all information in the operating instructions and complying with the measures described therein is also a part of the intended use.

### 4 Remaining risks

The device is built according to the state-of-the-art and the recognized safety regulations.

If used and handled improperly, there could be life-threatening danger to the user or third parties, or the device could be impaired or there could be other property damage.

The device is only to be used for its intended purpose and only when it is in safe working condition.

Malfunctions which could affect safety must be corrected immediately.

## 5 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen			
Model	ROTIXA 500 RS			
Type	4950, 4950-50		4950-70, 4950-80	
Mains voltage ( $\pm 10\%$ )	230-240 V ~	220 V ~	230-240 V ~	220 V ~
Mains frequency	50 Hz	60 Hz	50 Hz	60 Hz
Connected load	3800 VA			
Current consumption	16 A			
Cooling medium	R 404A			
Max. capacity	4 x 1000 ml			
Allowed density	1.2 kg/dm <sup>3</sup>			
Speed (RPM)	11500			
Force (RCF)	18038			
Kinetic energy	59620 Nm			
Obligatory inspection (BGR 500)	yes			
Ambient conditions (EN / IEC 61010-1)	<p>Indoors only Up to 2000 m above sea level 5°C to 35°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.</p>			
– Set-up site				
– Altitude				
– Ambient temperature				
– Humidity				
– Excess-voltage category (IEC 60364-4-443)	II			
– Pollution degree	2			
Device protection class	I			
Not suitable for use in explosion-endangered areas.				
EMC	EN / IEC 61326-1, Class B			
– Emitted interference, Interference immunity				
Noise level (dependent on rotor)	$\leq 65$ dB(A)		$\leq 63$ dB(A)	
Dimensions	<p>650 mm 814 mm 973 mm</p>			
– Width				
– Depth				
– Height	973 mm			
Weight	approx. 219 kg		approx. 233 kg	

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen	
Model	ROTIXA 500 RS	
Type	4950-08, 4950-58	4950-78, 4950-88
Mains voltage ( $\pm 10\%$ )	208 V ~	
Mains frequency	60 Hz	
Connected load	3800 VA	
Current consumption	18 A	
Cooling medium	R 404A	
Max. capacity	4 x 1000 ml	
Allowed density	1.2 kg/dm <sup>3</sup>	
Speed (RPM)	11500	
Force (RCF)	18038	
Kinetic energy	59620 Nm	
Obligatory inspection (BGR 500)	yes	
Ambient conditions (EN / IEC 61010-1)	<p>Indoors only Up to 2000 m above sea level 5°C to 35°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.</p>	
– Set-up site		
– Altitude		
– Ambient temperature		
– Humidity		
– Excess-voltage category (IEC 60364-4-443)	II	
– Pollution degree	2	
Device protection class	I	
Not suitable for use in explosion-endangered areas.		
EMC	FCC Class B	
– Emitted interference, Interference immunity		
Noise level (dependent on rotor)	$\leq 65$ dB(A)	$\leq 63$ dB(A)
Dimensions		
– Width	650 mm	
– Depth	814 mm	
– Height	973 mm	
Weight	approx. 225 kg	approx. 239 kg

## 6 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.


- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
- The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.
- With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.
- For centrifuges without temperature control, when the room temperature is increased and/or if the device is frequently used, the centrifuging chamber could be heated up. Therefore, it can't be ruled out that the sample material might be changed due to the temperature.

- Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). Before centrifuge vessels are used which are not listed in the chapter "Appendix, Rotors and accessories", the user must make sure they can be used by asking the manufacturer.
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
  - inflammable or explosive materials
  - materials that react with one another producing a lot of energy.
 When centrifuging petrochemical samples according to the specifications of the ASTM, the user must carry out a risk assessment.

- When centrifuging hazardous substances or mixtures, which are toxic, radioactive or contaminated with pathogenic microorganisms, suitable measures must be taken by the user.  
Fundamentally, centrifuge containers with special screw closures must be used for hazardous substances. For materials of risk groups 3 and 4, in addition to sealable centrifuge containers, a bio-safety system must be used (see the "Laboratory Bio-safety Manual" from the World Health Organization).  
In a bio-safety system, a bio-seal (sealing ring) prevents droplets and aerosols from escaping.  
If the hanger of a bio-safety system is used without the lid, the sealing ring must be removed from the hanger to prevent damage to the sealing ring during the centrifugation run.  
Damaged bio-safety systems are no longer microbiologically sealed.  
If a bio-safety system is not used, a centrifuge is not microbiologically sealed for the purposes of the standard EN / IEC 61010-2-020.  
When closing a bio-safety system, follow the instructions in the chapter "Closing biosafety systems".  
For the available bio-safety systems, see the chapter "Appendix, Rotors and accessories". If in doubt, you can get the information you need from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:  
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
  - the centrifuge is operated in accordance with the operating instructions,
  - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations,
  - the tests for device safety required in the respective countries, e.g. in Germany in acc. with DGUV (Deutsche Gesetzliche Unfallversicherung or Engl.: German Statutory Accident Insurance) regulation 3, are carried out by an expert.

## 7 Unpacking and setting up the centrifuge

- Remove the packaging.
- Remove the front wooden beam (a).
- Fasten the metal rails (b) to the wooden pallet with two nails each.
- Push the front wooden beam (a) under the metal rails (b) to support these.
- Place the open-end wrench included in delivery (size 10 mm) on the surfaces (d) and turn the feet of the device (c) upward.
- Roll the centrifuge carefully down from the wooden pallet using the metal rails (b).
- Push the centrifuge to its set-up location.
- Place the open-end wrench included in delivery (size 10 mm) on the surfaces (d) and turn the device feet (c) downward until the castor wheels (f) no longer have contact with the floor.
- Align the centrifuge horizontally by turning the device feet (c).
- Turn the hexagon nuts (e) upward with the open-end wrench included in delivery (size 19 mm) and tighten to secure the position of the device feet (c).

## 8 Delivery checklist

- 1 Power cord
- 1 Single open-end wrench, size 10 mm (for turning the device feet)
- 1 Double open-end wrench, size 17/19 mm (for tightening the hexagon nuts on the device feet)
- 1 Square wrench (for installing and removing the rotor)
- 1 Allen key (to unlock the lid in an emergency)
- 1 Lubricating grease for trunnions
- 1 Tube emptier
- 10 Cover caps ( $\varnothing$ 12 mm) for Torx screw heads, size T20
- 1 Operating manual
- 3 Program-data sheets

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

## 9 Transport and storage

### 9.1 Transport

When the device and accessories are transported, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 20% to 80%, non-condensing

### 9.2 Storage

The device and the accessories may only be stored in closed and dry rooms.

When the device and accessories are stored, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 20% to 80%, non-condensing

## 10 Initial operation

- Only for types 4950-08, 4950-58, 4950-78 and 4950-88:  
According to the laboratory device standard EN / IEC 61010-1, there must be a switch installed in the building installation to disconnect the power to the device.  
This switch has to be placed close to the centrifuge, must be easily reachable for the user and has to be characterized as a separator for this unit. The switch must also be able to be secured against being switched on again.

- Types 4950-70, 4950-78, 4950-80 and 4950-88 may only be installed by authorized skilled personnel.

These devices must be installed according to the enclosed set-up and installation instructions, AH4950EN.

- **Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

**When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.**

- Ventilation openings may not be blocked.  
A distance of 300 mm must be maintained from the ventilation slots and openings of the centrifuge.
- Centrifuge with PE connector:  
If required, connect the PE connector at the back of the device with an additional medical potential equalization system.
- Centrifuge with optical interface:  
Connect the optical centrifuge interface with a light-conducting cable (not included in the scope of delivery) to the PC.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the power cord to a standard mains socket. For connection ratings refer to Chapter "Technical specification".
- Turn on the mains switch. Switch position "I".  
The following display follows:  
1. Type of centrifuge, 2. The last rotor code detected by the rotor identification and the maximum speed of the rotor (n-max-Rotor), 3. The programme version, 4. **⚡ OPEN ⚡ OEFFNEN**.
- Open the lid.  
The centrifuge data of the last used program or program 1 will be displayed.

## 11 Interface (only for centrifuge with interface)

As an option, the device can be equipped with an optical interface.

The optical interface is labeled with the  symbol.

The centrifuge can be controlled and data queried via this interface.

## 12 Opening and closing the lid


### 12.1 Opening the lid

 The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".


- Press the button  on the front panel. The lid unlocks via the motor, the illuminated pushbutton  extinguishes and in the display the symbol  changes to the symbol .

### 12.2 Closing the lid

 Do not put your fingers between lid and housing.  
Do not bang the lid shut.

- Place the lid and slightly press it down at the tubular handle. The locking action is effected by motor. The button  lights up and the symbol  in the display turns to the symbol .

## 13 Installation and removal of the rotor


- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The motor shaft dog (D) has to fit in the rotor slot (B). The alignment of the groove is labelled on the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.

## 14 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.
- With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position. Hangers that are marked with a set number (e.g. S001/4) may only be used in the set.
- The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".  
In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Rotor is evenly loaded


**Not permitted!**

Rotor is not evenly loaded

- On certain hangers, the weight of the maximum load or the weight of the maximum load and the maximum weight of the completely loaded hanger is specified. This weight may not be exceeded. In case of exception, see chapter "Centrifugation of materials or mixtures of materials with a density higher than  $1.2 \text{ kg/dm}^3$ ". The weight specification of the maximum load includes the total weight of the adapter, centrifuge container and contents.
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- The centrifuge containers may only be filled outside of the centrifuge.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.

In the case of angle rotors, the centrifuging vessels may only be filled so far that no fluid can be expelled from them while the centrifuge is running.

- When loading the angle rotors, no liquid may get into the angle rotors or centrifuging chamber.
- When loading the hangers and swing-out rotors and when swivelling out the hangers while the centrifuge is running, no liquid may enter the hangers or the centrifuging chamber.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.

## 15 Closing biosafety systems


In order to guarantee that it is sealed, the lid to the biosafety system has to be firmly closed.

In order to prevent the packing ring from twisting when opening and closing the cover, the packing ring must be lightly rubbed with talcum powder or a rubber care product.


For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.

### Lid with screw top and hole in turning handle


- Place the lid centred on the rotor.
- Insert the supplied key into the hole in the handle and close the lid firmly by turning the key in a clockwise direction.

### Lid with screw closure

- Place the lid on the hanger.
- Tightly seal the lid by hand by turning it clockwise.

## 16 Control and display devices

### 16.1 Status symbols

Lid open.

Lid closed.

Rotation indicator. The indication appears while the centrifuge is running as long as the rotor is turning.

Centrifugation run stopped or finished. Indication after completion of the centrifugation run as long as the rotor is turning. The indicator flashes following an emergency stop.

**LOCK 1, LOCK 2** Keyswitch position.

**LOCK 4, LOCK 5** Program interlocking when serial communication provided (for centrifuges with serial communication only).

**PC, PE** Serial communication (for centrifuges with serial communication only).

Operating errors or faults that occur are shown symbolically in the display (see Chapter "Faults").

### 16.2 Control knob

For setting the individual parameters.  
Turning anticlockwise reduces the value. Turning clockwise increases the value.

### 16.3 Keys and set-up options

Running time, parameter t/min:sec  
 1. Parameter t/min: Adjustable from 1 - 999 min, in 1 min steps.  
 2. Parameter t/ :sec Adjustable from 1 - 59 s, in 1 second steps.  
 3. Continuous running "---:--"

Interrogation of the integral RCF, parameter JRCF.

1. Speed, parameter RPM  
 A numerical value from 50 RPM up to the maximum speed of the rotor (n-max-Rotor) can be set. Maximum speed of the rotor, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".  
 Adjustable in steps of 10.  
 2. Interrogation of the maximum speed of the rotor, parameter n-max-Rotor.

1. Relative centrifugal force, parameter RCF/RZB  
 A numerical value can be set, which gives a speed between 50 RPM and the maximum speed of the rotor (n-max-Rotor). Adjustable in steps of 1.  
 2. Interrogation of the maximum RCF of the rotor, parameter RCF-max-Rotor.

Run-up parameters  
 1. Run-up steps, parameter  $\curvearrowright_{1-9}$ . Step 9 = shortest run-up time, ... Step 1 = longest run-up time.  
 2. Run-up time, parameter  $\curvearrowright_{\text{min:sec}}$ . The time range that can be set is dependent upon the set speed.

Run-down parameters  
 1. Braking steps, parameter  $\curvearrowleft_{0-9}$ . R = linear braking curve, B = similar to an exponential braking curve.  
 Step R9, B9 = short run-down time, ... Step R1, B1 = long run-down time, Step R0 = unbraked run-down.  
 2. Run-down time, parameter  $\curvearrowleft_{\text{min:sec}}$ . The time range that can be set is dependent upon the set speed.  
 3. Brake de-energisation speed, parameter  $n^{(*)}$ /RPM  
 After reaching this speed, unbraked run-down occurs.

1. Temperature setpoint, parameter T/°C. Adjustable from -20°C to +40°C, in 1°C steps (heating/cooling option settable between -20°C and +60°C). The lowest temperature that can be achieved is dependent upon the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").  
 2. Centrifuging radius, parameter r/mm. Entry in mm. Centrifuging radius, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".


- START 1. Start centrifugation run. The rotation indicator  appears.  
2. Acceptance of changes during the centrifugation run.
- STOP Stop centrifugation run.  
The rotor decelerates with the preset run-down parameters. Pressing the key twice will initiate an emergency stop.
- PROG Select program location, parameter PROG-Nr.
- RCL Calling up of programs.
- STO Saving of programs.  
89 programs can be saved (program locations 1 to 89). Note: The program locations "----" and 90 to 99 are used as automatic buffers (see Chapter "Automatic buffer"). No programs can be saved in these program locations.

## 17 Entering centrifugation parameters

 It is only possible to enter a parameter when the input field is displayed inversely (dark background). An inverse input field will go off automatically after 10 seconds.

### 17.1 Running time


#### 17.1.1 Time selection

- Press the  key repeatedly until the parameter t/min: or t/ :sec is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

#### 17.1.2 Continuous running

- Select the parameters t/min: and t/ :sec one after the other (see Chapter "Time selection"), and set both to zero using the control knob . "---:--" will be displayed in the input field.

### 17.2 Speed (RPM)

- Press the  key repeatedly until the parameter RPM is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

### 17.3 Relative centrifugal force (RCF/RZB)

- Press the  key repeatedly until the parameter RCF/RZB is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

The centrifuging radius indicator will flash while the adjustment is being made.

### 17.4 Run-up and run-down parameters


The run-up and run-down parameters that have been set are displayed.

x: 1-9 = run-up step, t = run-up time

y: R1-R9, B1-B9 = braking step, R0 = unbraked run-down, t = run-down time, n<sup>(\*)</sup> = brake de-energisation speed

#### 17.4.1 Run-up step

- Press the  key repeatedly until the parameter <sub>1-9</sub> or  is displayed. The input field will be displayed inversely.
- Set the desired step by means of the control knob .

#### 17.4.2 Run-up time

- Press the  key repeatedly until the parameter <sup>min:sec</sup> is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

If a run-up time is set that is longer than the running time, the centrifugation run will end before the set speed is reached.

### 17.4.3 Braking step

- Press the  $\square$  key repeatedly until the parameter  $\backslash_{0-9}$  or  $\backslash$  is displayed. The input field will be displayed inversely.
- Set the desired step by means of the control knob  $\odot$ .

B-braking steps can only be set for special rotors.

### 17.4.4 Run-down time

If a brake de-energisation speed has been set, it is not possible to set a run-down time.

- Press the  $\square$  key repeatedly until the parameter  $\backslash_{\text{min:sec}}$  is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob  $\odot$ .

### 17.4.5 Brake de-energisation speed

- Press the  $\square$  key repeatedly until the parameter  $n^{(+)}/\text{RPM}$  is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob  $\odot$ .

## 17.5 Radius/Temperature

### 17.5.1 Radius

- Press the  $\square$  key repeatedly until the parameter  $r/\text{mm}$  is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob  $\odot$ .

The value of the RCF/RZB adapts automatically when the radius is changed; this is shown by the indication flashing.

### 17.5.2 Temperature

- Press the  $\square$  key repeatedly until the parameter  $T/^\circ\text{C}$  is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob  $\odot$ .

### 17.6 Automatic buffer

The buffer comprises the program locations "----" and 90 to 99. Every time a centrifugation run is started, modified centrifugation data are automatically stored in program location "----". The modified centrifugation data for the last 11 centrifugation runs are stored in the buffer and can be called up (see Chapter "Calling up programs").

## 18 Programming

### 18.1 Entering / modifying programs

- Set the desired parameters (see Chapter "Entering centrifugation parameters").
- Press the  $\text{PROG}$  key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
- Set the desired program location by means of the control knob  $\odot$ .  
If the program location indicator flashes, then this program location has already had centrifugation data assigned to it. In this case, set a free program location or overwrite the centrifugation data by continuing.
- Press the  $\text{STO}$  key once to save the settings in the desired program location.  
Press the  $\text{STO}$  key twice to overwrite centrifugation data already stored.

### 18.2 Calling up programs

- Press the  $\text{PROG}$  key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
- Set the desired program location by means of the control knob  $\odot$ .
- Press the  $\text{RCL}$  key. The centrifugation data for the selected program location will be displayed.

## 19 Centrifugation

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, and **IMBALANCE / UNWUCHT** is displayed.

If the speed in the selected program is higher than the maximum speed of the rotor (n-max-Rotor), it will not be possible to start a centrifugation run. The message **N > ROTOR MAX 96** will be displayed (see Chapter "Faults").

A centrifugation run can be stopped at any time by pushing the key **STOP**.

It is possible to switch over between RPM and RCF display at any time using the **RPM** and **RCF** keys. When working with the RCF display, it is necessary to enter the centrifuging radius.

If **OPEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

If **R xx n-max xxxxx** is displayed, then no centrifugation run has taken place as the rotor was changed beforehand, refer to Chapter "Rotor Identification".




- Switch on the mains supply switch. Switch position I.
- Load the rotor and close the centrifuge lid.

### 19.1 Centrifugation with time selection

- Set the time or call up a program with a pre-set time (see Chapter "Entering centrifugation parameters" or "Calling up programs").
- Press the **START** key. The rotation indicator  will appear as long as the rotor is turning.
- When the time has expired or if the centrifugation run is interrupted by pressing the **STOP** key, run-down will take place with the selected run-down parameters. When the rotor comes to a stop, the  symbol will flash in the display until the lid is opened.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature and the remaining time will be displayed.

### 19.2 Continuous running

- Select the  symbol or call up a continuous running program (see Chapter "Entering centrifugation parameters" or "Calling up programs").
- Press the **START** key. The rotation indicator  will appear as long as the rotor is turning. The time count starts from 00:00.
- Press the **STOP** key to stop the centrifugation run. Run-down will take place with the selected run-down parameters. When the rotor comes to a stop, the  symbol will flash in the display until the lid is opened.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature and the time run will be displayed.

## 20 Changing settings during the centrifugation run.

The running time, the speed, the relative centrifugal force (RCF/RZB), the run-up and run-down parameters as well as the temperature can be changed while the centrifuge is running.

The parameters can only be changed one at a time and in succession.

- Select the desired parameter and change the value using the control knob  (see Chapter "Entering centrifugation parameters").
- Press the **START** key. The modified setting will be stored in program location "----" (see Chapter "Automatic buffer"). The original program will not be overwritten.

## 21 Integral RCF (∫RCF)

The integral RCF (∫RCF) is a measure for the sedimentation effect ( $\int n^2 dt$ ). This numerical value is used to compare centrifugation runs. To interrogate the integral RCF, press and hold the **RCF** key.

## 22 Displaying the maximum speed of the rotor

- Press the **RPM** key repeatedly until the parameter RPM is displayed and the input field is shown inversely.
- Press the **RPM** key once more and hold it pressed; the maximum speed of the rotor (n-max-Rotor) will be displayed.

### 23 Displaying the maximum RCF of the rotor

- Press the **[RCF]** key repeatedly until the parameter RCF/RZB is displayed and the input field is shown inversely.
- Press the **[RCF]** key once more and hold it pressed; the maximum RCF of the rotor (RCF-max-Rotor) will be displayed.

### 24 Emergency stop

- Press the **[STOP]** key twice. The **STOP** symbol will flash in the display.

In the case of an emergency stop, run-down takes place with braking step R9 (shortest run-down time). If braking step R0 was selected, the run-down time will be longer than with braking step R9 for technical reasons.

### 25 Audible signal

The audible signal sounds in accordance with the following philosophy:

<b>OFF</b>	<ul style="list-style-type: none"> <li>• at 2 s intervals on the occurrence of a fault</li> </ul>
<b>ON1</b>	<ul style="list-style-type: none"> <li>• at 2 s intervals on the occurrence of a fault</li> <li>• at 30 s intervals on completion of a centrifugation run and the rotor coming to rest.</li> </ul>
<b>ON2</b>	<ul style="list-style-type: none"> <li>• at 2 s intervals on the occurrence of a fault</li> <li>• at 30 s intervals on completion of a centrifugation run and the rotor coming to rest.</li> <li>• every time a key is pressed</li> </ul>

The audible signal is silenced by opening the lid or pressing any key.

The signal can be set as follows when the rotor is at standstill.

- Open the lid.
- Press and hold the **[ ]** key for 8 seconds.  
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Set the desired function (OFF, ON1, ON2) by means of the control knob **○**.
- Press the **[START]** key to save the setting.  
**\*\*\* ok \*\*\*** will be displayed briefly as a confirmation.

### 26 Interrogating the operating hours

It is only possible to interrogate the operating hours when the rotor is at standstill.

- Open the lid.
- Press and hold the **[ ]** key for 8 seconds.  
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Press the **[ ]** key once more.  
The operating hours (**CONTROL:** ) for the centrifuge will be displayed.  
The operating hours display will go off automatically after 10 seconds.

### 27 Setting the date and time

It is only possible to set the date and time when the rotor is at standstill.

- Open the lid.
- Press and hold the **[ ]** key for 8 seconds.  
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Press the **[ ]** key twice more.  
The date and time will be displayed (a: year, mon: month, d: day, h: hours, min: minutes).
- Press the **[ ]** key repeatedly until the input field for the desired parameter is shown inversely.
- Set the desired value by means of the control knob **○**.
- Press the **[START]** key to save the setting.  
**\*\*\* ok \*\*\*** will be displayed briefly as a confirmation.
- To exit the date and time set-up, press any key apart from the **[ ]**, **[ ]** and **[START]** keys.

### 28 Centrifugation data displayed after switching on.

After switching on, the centrifugation data for program 1 or for the program last used will be displayed. This can be set as follows:

- Open the lid.
- Switch the mains supply switch off and on again. Switch position I.
- Press the **[STOP]** key when the first visual change appears in the display (inverse display).  
**PROGRAM 1, LAST PROGRAM** will be displayed.
- Set the desired function by means of the control knob **○**.
- Press the **[START]** key to save the setting.  
**\*\*\* ok \*\*\*** will be displayed briefly as a confirmation.

## 29 Immediate display of the centrifugation data after switching on

- Switch on the mains supply switch. Switch position I.
- Press any key apart from the **STOP** key when the first visual change appears in the display (inverse display). The centrifugation data will be displayed immediately.

## 30 Key-operated switch.

The following program interlocks can be set by the key-operated switch:

Left key position: **LOCK 1** will be displayed.

Programs can be called up but not changed.

Right key position: **LOCK 2** will be displayed.

No programs can be called up or changed.

Middle key position: no status display.

No program interlock. Programs can be called up and changed.

## 31 Program linking (for centrifuges with program linking only).

Several centrifugation runs can be linked together with the help of the program linking function.

### 31.1 Linking programs

Program linking is only possible with programs in which run-up and braking steps have been set.

Before linking, the programs must be stored in the desired sequence either by entering the program or by calling up the program (see Chapter "Programming").

The program locations must be consecutive (e.g. program locations 10+11+12).


1. Press the **PROG** key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
2. Set the program location for the first program (XX+) by means of the control knob  $\odot$ .
3. Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
4. Press the **PROG** key twice in order to select the parameter PR-PART. The input field will be displayed inversely.
5. Press the **STO** key twice. The program will be linked and the program number for the next program location (+XX+) will be displayed.
6. Press the **RCL** key twice. The centrifugation data for the selected program location will be displayed.
7. Press the **STO** key twice. The program will be linked and the program number for the next program location (+XX+) will be displayed.
8. Keep repeating steps 6 and 7 until all programs have been linked.
9. Press the **PROG** key to exit. The program number for the last program (+XX) will be displayed.

### 31.2 Changing program links

- Call up the desired program (see Chapter "Calling up programs"), change the desired parameters (see Chapter "Entering centrifugation parameters"), and save the modified centrifugation data to the same program location (see Chapter "Entering / modifying programs"). As a result of saving, the program linking will be removed.
- Re-link the programs (see Chapter "Linking programs").

### 31.3 Centrifugation run with program linking

- Press the **PROG** key twice in order to select the parameter PR-PART. The input field will be displayed inversely.
- Set the program location for the first program (XX+) by means of the control knob  $\odot$ .
- Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
- Press the **START** key. The rotation indicator  $\odot$  will appear as long as the rotor is turning.  
The run-up and braking steps for the program linking will be displayed:

- When the time in the last program has expired, run-down will take place with the braking step of the last program. If the centrifugation run is interrupted by pressing the **STOP** key, run-down will take place with the braking step of the program running at the time.

### 31.4 Deleting program linking

- Press the **PROG** key in order to display the parameter PROG-Nr. The input field will be displayed inversely.
- Set the program location for the first program (XX+) by means of the control knob  $\odot$ .
- Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
- Press the **PROG** key twice in order to display the parameter PR-PART. The input field will be displayed inversely.
- Press the **STO** key twice.
- Press the **PROG** key.

## 32 Cooling

The temperature set-point can be adjusted from -20°C to +40°C. On centrifuges with a heating/cooling option, the temperature set-point is settable from -20°C to +60°C. The lowest obtainable temperature is dependent on the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). If the actual temperature and the set temperature deviate more than 5 °C, the temperature value in the display will start to flash on and off.

### 32.1 Standby-cooling

With rotor standstill and closed lid the centrifugal chamber is cooled to the pre-selected temperature. The temperature set-point is shown in the display.

### 32.2 Pre-cooling the rotor

For rapid pre-cooling of the unloaded rotor and the accessories, it is recommended that the centrifuge is run with the continuous running settings and a speed of ca. 20% of the maximum rotor speed.

## 33 Heating

During operation of the centrifuge the centrifuge chamber can be heated as required to the preselected temperature. The heating is switched off when the rotor is stationary.

Risk of burning! The surface temperature of the heating element in the centrifugal chamber of the centrifuge can reach 500°C / 932°F. Do not touch the heating element.

### 34 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left( \frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge container. For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

### 35 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm<sup>3</sup>

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm<sup>3</sup>

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

In the exceptional case that the maximum loading indicated on the hanger is exceeded, the speed must also be reduced.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, maximum load 300 g, actual load 350 g

$$n_{\text{red}} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} \times 4000 \text{ RPM} = 3703 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

### 36 Rotor recognition

Rotor recognition is carried out after each start of the centrifugation run.

If the rotor has been changed, the centrifugation run is canceled after rotor recognition. The rotor code (R xx) and the maximum speed (n-max=xxxxx) of the rotor are displayed.

A further operation of the centrifuge is only possible after a single opening of the lid.

If the maximum speed of the rotor being used is less than the set speed, the speed is limited to the rotor's maximum speed.


### 37 Emergency release

During a power failure the lid cannot be unlocked by motor. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains. In the case of permanently connected devices, switch off the mains switch to disconnect the mains supply of the device in the building installation and secure it against being switched on again, e.g. by locking the switch.

Open the lid only during rotor standstill.

- Switch off the mains switch (switch setting "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Insert the hexagon socket head wrench horizontally into the drilling (A) and carefully rotate by half a turn in clockwise direction until the lid can be opened.
- Pull the hexagon socket head wrench out of the drilling again.
- After turning the centrifuge on again, press the button  so that the motor-driven lid locking once again assumes the normal position (opened).

### 38 Maintenance and servicing

The device can be contaminated.

Pull the mains plug before cleaning. In the case of permanently connected devices, switch off the mains switch to disconnect the mains supply of the device in the building installation and secure it against being switched on again, e.g. by locking the switch.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
  - have a pH between 5 - 8
  - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

#### 38.1 Centrifuge (housing, lid and centrifuging chamber)

##### 38.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

### 38.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants:  
ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

### 38.1.3 Removal of radioactive contaminants


- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:  
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

## 38.2 Rotors and Attachments

### 38.2.1 Cleaning and care

- In order to avoid corrosion and changes in materials, the rotors and accessories have to be cleaned regularly with soap or with a mild cleaning agent and a moist cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents:  
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- The sealing rings of bio-safety systems must be cleaned weekly.  
The sealing rings are made of silicone. To guarantee the leak-tightness of the bio-safety systems, the sealing rings must not be handled with talcum powder after cleaning or autoclaving.  
Each time before using the bio-safety system, all parts of the bio-safety system must be visually inspected for damage. In addition, the sealing ring(s) of the bio-safety system must be checked to make sure they are in the correct installation position.  
Damaged parts of the bio-safety system must be exchanged immediately.  
If there are signs of crack formation, brittleness or wear, immediately replace the sealing ring in question. In the case of lids with sealing rings which cannot be replaced, the entire lid must be exchanged.  
For the available bio-safety systems, see the chapter "Appendix, Rotors and accessories".
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and accessories have to be checked weekly for wear and corrosion. For swing-out rotors, it is important to check the area of the lifting lugs, for hangers, the grooves and the base should be checked for cracks.

Example: Crack in the groove area:

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

### 38.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants:  
ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

### 38.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:  
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

### 38.2.4 Trunnions

With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.

### 38.2.5 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

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For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

### 38.2.6 Calculation of the number of run cycles carried out

To calculate the number of carried out run cycles (centrifugation runs), the run-time per centrifugation run and the operating hours of the device are required.

To query the operating hours, see the chapter "Operating hours query".

If centrifugation runs were carried out with different run-times, the shortest run-time is to be used for the calculation.

The number of carried out run cycles (centrifugation runs) are calculated as follows:

Number of executed running cycles = operating hours [h] × 60 ÷ run-time [min]

e.g.: Operating hours 2000 h, run-time 5 min

Number of executed running cycles = 2000 × 60 ÷ 5 = 24000

### 38.3 Autoclaving

The following accessories may be autoclaved at 121°C / 250°F (20 min):

- Swing-out rotors
- Angle rotors made of aluminium
- Hangers made of metal
- Lid with bio-seal
- Adapter

Nothing definitive can be said about the degree of sterility.

The lids of the rotors and containers must be removed before autoclaving.

Autoclaving accelerates the ageing process of plastics. In addition, it can cause discolourations in plastics.

After autoclaving, the rotors and the accessories must be visually inspected for damage and any damaged parts must be exchanged immediately.

If there are signs of crack formation, brittleness or wear, immediately replace the sealing ring in question.

In the case of lids with sealing rings which cannot be replaced, the entire lid must be exchanged.

To guarantee the leak-tightness of the bio-safety systems, the sealing rings must not be handled with talcum powder after autoclaving.

### 38.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts, glass splinters and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splinters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

### 39 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

Perform a MAINS RESET:

- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "I").

Message / fault		Cause	Remedy
No display	---	No voltage. Overvoltage protection tripped out.	<ul style="list-style-type: none"> <li>- Check supply voltage.</li> <li>- Mains switch ON.</li> </ul>
TACHO - ERROR	01	Faulty speedometer.	<ul style="list-style-type: none"> <li>- Open the cover.</li> <li>- Switch off the mains switch (switch position "0").</li> <li>- Wait at least 10 seconds.</li> <li>- Turn the rotor vigorously by hand.</li> <li>- Switch on the mains switch again (switch position "I"). The rotor must turn during switch-on.</li> </ul>
	02	No rotor installed. Defective motor, frequency converter or drive.	
IMBALANCE / UNWUCHT	---	The rotor is unevenly loaded.	<ul style="list-style-type: none"> <li>- Open lid.</li> <li>- Check the loading of the rotor, see chapter "Loading the rotor".</li> <li>- Repeat the centrifugation run.</li> </ul>
CONTROL - ERROR	04, 06 - 09	Error in lid locking or lid closure.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET.</li> </ul>
N > MAX	05	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10	Incorrect rotor coding	
MAINS INTERRUPT	---	Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> <li>- Open lid.</li> <li>- Push <b>START</b> button.</li> <li>- Repeat the centrifugation run if necessary.</li> </ul>
VERSIONS-ERROR	12	Mismatch between electronic components	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET.</li> </ul>
SER I/O - ERROR	30 - 38	Error / defective interface	
° C * - ERROR	50 - 56, 58	Error / defective cooling	
LOCK - ERROR	57	Error / defective program locking	
FU / CCI - ERROR	60 - 83	Error / defective motor control	
CONTROL-ERROR	26, 90 - 95, 97 - 99	Error / defective control unit	
N > ROTOR-MAX	96	Speed in the selected program greater than the maximum speed of the rotor (n-max-Rotor).	<ul style="list-style-type: none"> <li>- Check the set speed.</li> <li>- Reduce the set speed</li> </ul>
		The rotor has been changed. The maximum speed of the installed rotor is higher than that of the rotor that was previously used and this has not yet been detected by rotor recognition.	<ul style="list-style-type: none"> <li>- Set a speed up to the maximum speed of the rotor that was previously used. Press the <b>START</b> key to run rotor recognition. See chapter, "Rotor Recognition".</li> </ul>

#### 40 Returning Devices

If the device or its accessories are returned to Andreas Hettich GmbH & Co. KG, in order to provide protection for people, the environment and materials, it has to be decontaminated and cleaned before being shipped.

We reserve the right to refuse contaminated devices or accessories.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

#### 41 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)  
 3) not possible to close the lid 5053  
 10) with decanting aid













4296	5051 + 5053									
<p><b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b></p>  <p>∠ 90°</p>	 									
	max. Laufzyklen / max. cycles: 50000 max.Beladung / max. load: 500 g									
	5259	6306	5248	5264		5267		5281		
	 <sup>3)</sup>	 <sup>3)</sup>								
0513	0509				0500			2078	0536	
										
Kapazität / capacity ml	50	15	8,5 - 10	4 - 7	9	1,1-1,4	3	1,5	2,0	
Maße / dimensions Ø x L mm	29 x 115	17 x 120	16 x 100	16 x 75	14 x 100	8 x 66	10 x 60	11 x 38		
Anzahl p. Rotor / number p. rotor	12	42	72	72	72	120	120	96		
Drehzahl / speed RPM	4000	4000	4000	4000	4000	4000	4000	4000		
RZB / RCF	3363	3434	3291	3309	3309	3274	3274	3363		
Radius / radius mm	188	192	184	185	185	183	183	188		
 9 (97%) sec					33					
 9 sec					50					
Temperatur / temperature °C <sup>1)</sup>	0									

4296	5051 + 5053									
<p><b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b></p>  <p>∠ 90°</p>	 									
	max. Laufzyklen / max. cycles: 50000 max.Beladung / max. load: 500 g									
	5268									
										
										
Kapazität / capacity ml	1,6 - 5	4 - 7	5	2,6 - 3,4	4,9					
Maße / dimensions Ø x L mm	13 x 75	13 x 100	13 x 75	13 x 65	13 x 90					
Anzahl p. Rotor / number p. rotor	72	72	72	72	72					
Drehzahl / speed RPM	4000	4000	4000	4000	4000					
RZB / RCF	3345	3345	3345	3345	3345					
Radius / radius mm	187	187	187	187	187					
 9 (97%) sec					33					
 9 sec					50					
Temperatur / temperature °C <sup>1)</sup>	0									

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







1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)  
 3) not possible to close the lid 5053




4296		5051 + 5280						5053	
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>  									
		1662						1670	
								 4)	
		1663	1664	1665	1666	1667	1668	1663	1664
									
Kapazität / capacity	ml	1	2	4	8	3 x 2	4 x 1	1	2
Maße / dimensions $\varnothing / A$	mm <sup>2</sup>	6,2 / 30	8,7 / 60	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30	6,2 / 30	8,7 / 60
Anzahl p. Rotor / number p. rotor		12	12	12	12	12	12	12	12
Filterkarten / filter cards		1675	1675	1675	1676	1677	1668	1692	1692
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		2290 / 3274	2290 / 3274	2290 / 3274	2290 / 3274	2290 / 3274	2290 / 3274	2290 / 3274	2290 / 3274
Radius / radius	mm	128 / 183	128 / 183	128 / 183	128 / 183	128 / 183	128 / 183	128 / 183	128 / 183
 9 (97%)	sec							33	
 9	sec							50	
Temperatur / temperature	°C <sup>1)</sup>							0	

4296		5051 + 5280						5053	
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>  									
		1670				1470			
		 4)							
		1665	1666	1667	1668	1475	1471		
									
Kapazität / capacity	ml	4	8	3 x 2	4 x 1	2 x 8	1 x 8		
Maße / dimensions $\varnothing / A$	mm <sup>2</sup>	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30	17,5 / 240	17,5 / 240		
Anzahl p. Rotor / number p. rotor		12	12	12	12	12	12		
Filterkarten / filter cards		1692	1691	1694	1693	---	---		
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000		
RZB / RCF		2290 / 3274	2290 / 3274	2290 / 3274	2290 / 3274	2200 / 3202	2200 / 3202		
Radius / radius	mm	128 / 183	128 / 183	128 / 183	128 / 183	123 / 179	123 / 179		
 9 (97%)	sec							33	
 9	sec							50	
Temperatur / temperature	°C <sup>1)</sup>							0	

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


1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)  
 4) Object slide will not stand RCF values exceeding 1100



4296	5092 + 5093											
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    	 mit Bioabdichtung / with bio-containment 5)											
	5126		5125		5123		5124		5122		5128	
												
	0523		0526		0513		0521		0519		0553 0501	
												
Kapazität / capacity	ml	100	100	50	50	25	30	4	5	6		
Maße / dimensions	Ø x L	mm	40 x 115	44 x 100	29 x 115	34 x 100	24 x 100	25 x 110	12 x 60	12 x 75	12 x 82	
Anzahl p. Rotor / number p. rotor		6	6	12	6	24	24	72		72		
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000		4000		
RZB / RCF		3488	3488	3631	3488	3434	3434	3542		3542		
Radius / radius	mm	195	195	203	195	192	192	198		198		
 9 (97%)	sec	33										
 9	sec	50										
Temperatur / temperature	°C 1)	- 2										

4296	5092 + 5093											
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    	 mit Bioabdichtung / with bio-containment 5)											
	max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 500 g											
	5120					5121					5121-93	
												
0578								3)		0518		
												
Kapazität / capacity	ml	7	4,5 - 5	4 - 7	2,6 - 3,4	9 - 10	10	8,5 - 10	8	15	4 - 7	
Maße / dimensions	Ø x L	mm	12 x 100	11 x 92	13 x 100	13 x 65	16 x 92	15 x 102	16 x 100	16 x 125	17 x 100	16 x 75
Anzahl p. Rotor / number p. rotor		72	72	72	42	42	42	42	42	42	42	
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	4000	4000	
RZB / RCF		3542	3542	3542	3542	3542	3542	3542	3542	3542	3542	
Radius / radius	mm	198	198	198	198	198	198	198	198	198	198	
 9 (97%)	sec	33										
 9	sec	50										
Temperatur / temperature	°C 1)	- 2										

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











- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 3) not possible to close the lid 5053
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".






<p><b>4296</b></p> <p><b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b></p>  <p>∠ 90°</p>	<b>5092 + 5093</b>									
	mit Bioabdichtung / with bio-containment <sup>5)</sup>									
	max. Laufzyklen /max. cycles: 30000 max.Beladung / max. load: 500 g									
	5136									
										
		0518	2079							
Kapazität / capacity	ml	10	15	10	4 – 4,5	7,5 - 8,5	9 - 10	10	4 - 7	8,5 - 10
Maße / dimensions	∅ x L mm	16 x 80	17 x 100	17 x 70	15 x 75	15 x 92	16 x 92	15 x 102	16 x 75	16 x 100
Anzahl p. Rotor / number p. rotor		48	48	48	48	48	48	48	48	48
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		3488	3488	3488	3488	3488	3488	3488	3488	3488
Radius / radius	mm	195	195	195	195	195	195	195	195	195
9 (97%)	sec	33								
9	sec	50								
Temperatur / temperature	°C <sup>1)</sup>	- 2								

<p><b>4296</b></p> <p><b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b></p>  <p>∠ 90°</p>	<b>5092 + 5093</b>									
	mit Bioabdichtung / with bio-containment <sup>5)</sup>									
	max. Laufzyklen /max. cycles: 30000 max.Beladung / max. load: 500 g									
	5137									
										
Kapazität / capacity	ml	1,6 - 5	4 - 7	4,9	1,1 - 1,4	2,6 – 3,4	2,7 - 3	4,5 - 5	5	5
Maße / dimensions	∅ x L mm	13 x 75	13 x 100	13 x 90	8 x 66	13 x 65	11 x 66	11 x 92	13 x 75	13 x 75
Anzahl p. Rotor / number p. rotor		48	48	48	48	48	48	48	48	48
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		3488	3488	3488	3488	3488	3488	3488	3488	3488
Radius / radius	mm	195	195	195	195	195	195	195	195	195
9 (97%)	sec	33								
9	sec	50								
Temperatur / temperature	°C <sup>1)</sup>	- 2								

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







- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".





4296		5092 + 5093							5092			
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    												
		mit Bioabdichtung / with bio-containment <sup>5)</sup> max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 500 g										
		1791		5134	5135	5129		5138		6319	6319	
												
0530				0509				5127	---			
										 <sup>12)</sup>		
Kapazität / capacity	ml	250	25	50	15	1.1 – 1.4	2.7 – 3	2.6 – 3.4	1.6 – 5	250	290	
Maße / dimensions	∅ x L	mm	65 x 115	25 x 90	29 x 115	17 x 120	8 x 66	11 x 66	13 x 65	13 x 75	62 x 122	62 x 137
Anzahl p. Rotor / number p. rotor		6	18	12	42	72	72	72	72	6	6	
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	4000	4000	
RZB / RCF		3631	3363	3560	3631	3077	3077	3077	3077	3631	3631	
Radius / radius	mm	203	188	199	203	172	172	172	172	203	203	
 9 (97%)	sec							33				
 9	sec							50				
Temperatur / temperature	°C <sup>1)</sup>							- 2				

4294		4290 + 4291								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    										
		max. Laufzyklen / max. cycles: 30000 (4500 RPM); 45000 (4000 RPM), 60000 (3500 RPM) max. Beladung / max. load: 1200 g								
		mit Bioabdichtung / with bio-containment <sup>5)</sup>								
		4273  								
		0553	0501	0578						
										
Kapazität / capacity	ml	5	6	7	2,6 - 3,4	4,9	1,6 - 5	4 - 7		
Maße / dimensions	∅ x L	mm	12 x 75	12 x 82	12 x 100	13 x 65	13 x 90	13 x 75	13 x 100	
Anzahl p. Rotor / number p. rotor		200	200	200	200	200	200	200	200	
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	4551								
Radius / radius	mm	201								
 9 (97%)	sec	115								
 9	sec	116								
Temperatur / temperature	°C <sup>1)</sup>	2								

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






1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)  
 2) Centrifuge containers of glass will not stand RCF values exceeding 4000  
 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".










4294		4290 + 4291							
<p><b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b></p>  <p>↙ 90°</p>									
		max. Laufzyklen / max. cycles: 30000 (4500 RPM); 45000 (4000 RPM), 60000 (3500 RPM) max. Beladung / max. load: 1200 g							
		mit Bioabdichtung / with bio-containment 5)							
		4310	4311	4313	4314	4321	4338		
									
						0509	0513		
Kapazität / capacity	ml	12	10	9 - 10	50	15	50	10	
Maße / dimensions	∅ x L mm	17 x 100	15 x 102	16 x 92	29 x 115	17 x 120	29 x 115	16 x 80	
Anzahl p. Rotor / number p. rotor		112	132	132	32	68	32	168	
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	
RZB / RCF	2)	4437	4573	4573	4528	4618	4618	4551	
Radius / radius	mm	196	202	202	200	204	204	201	
	9 (97%)						115		
	9						116		
Temperatur / temperature	°C 1)						2		

4294		4290 + 4291							
<p><b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b></p>  <p>↙ 90°</p>									
		max. Laufzyklen / max. cycles: 30000 (4500 RPM); 45000 (4000 RPM), 60000 (3500 RPM) max. Beladung / max. load: 1200 g							
		mit Bioabdichtung / with bio-containment 5)							
		4338							
									
		0500		0518					
									
Kapazität / capacity	ml	9	14	15	4 - 5,5	7,5 - 8,5	4 - 7	8,5 - 10	
Maße / dimensions	∅ x L mm	14 x 100	16 x 101	17 x 100	15 x 75	15 x 92	16 x 75	16 x 100	
Anzahl p. Rotor / number p. rotor		168	168	168	168	168	168	168	
Drehzahl / speed	RPM	4500							
RZB / RCF	2)	4551							
Radius / radius	mm	201							
	9 (97%)						115		
	9						116		
Temperatur / temperature	°C 1)						2		

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




- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".








4294		4290						
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    		  max. Laufzyklen /max. cycles: 30000 (4500 RPM); 45000 (4000 RPM), 60000 (3500 RPM) max.Beladung / max. load: 1200 g						
		4339	4323	4320				
								
								
Kapazität / capacity	ml	50	50	15				
Maße / dimensions	∅ x L	mm	29 x 115	29 x 115	17 x 120			
Anzahl p. Rotor / number p. rotor		40	40	112				
Drehzahl / speed	RPM	4500	4500	4500				
RZB / RCF	<sup>2)</sup>	4528	4618	4618				
Radius / radius	mm	200	204	204				
 9 (97%)	sec							115
 9	sec							116
Temperatur / temperature	°C <sup>1)</sup>							2

4294		4295-A + 4229-B								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    		   max. Laufzyklen /max. cycles: 15000 (4500 RPM); 50000 (4000 RPM) max.Beladung / max. load: 1060 g								
		4226	4225	4224	4241	4245	4213			
										
			 							
Kapazität / capacity	ml	0.8	1.5	2.0	4	25	50	6	7	4.5 - 5
Maße / dimensions	∅ x L	mm	8 x 45	11 x 38	10 x 88	25 x 90	29 x 115	12 x 82	12 x 100	11 x 92
Anzahl p. Rotor / number p. rotor		312	336	252	32	32	32	192	192	192
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	<sup>2)</sup>	4777	3690 / 4887	4777	4777	4777	4958	4777	4777	4777
Radius / radius	mm	211	163 / 215	211	211	211	219	211	211	211
 9 (97%)	sec							115		
 9	sec							116		
Temperatur / temperature	°C <sup>1)</sup>							6		

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





1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)  
 2) Centrifuge containers of glass will not stand RCF values exceeding 4000








4294		4295-A + 4229-B							
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  ↙ 90°									
		max. Laufzyklen / max. cycles: 15000 (4500 RPM); 50000 (4000 RPM) max. Beladung / max. load: 1060 g							
		4213-93				4214		4214-93	
									
		0553	0501		0518				
									
Kapazität / capacity	ml	5	6	2,7 - 3	15	7,5 - 8,2	8,5 - 10	10	4 - 7
Maße / dimensions	∅ x L mm	12 x 75	12 x 82	11 x 66	17 x 100	15 x 92	16 x 100	15 x 102	16 x 75
Anzahl p. Gestell / number p. frame		48	48	48	30	30	30	30	30
Anzahl p. Rotor / number p. rotor		192	192	192	120	120	120	120	120
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	<sup>2)</sup>	4777	4777	4777	4777	4777	4777	4777	4777
Radius / radius	mm	211	211	211	211	211	211	211	211
 9 (97%)	sec	115							
 9	sec	116							
Temperatur / temperature	°C <sup>1)</sup>	6							

4294		4295-A + 4229-B							
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  ↙ 90°									
		max. Laufzyklen / max. cycles: 15000 (4500 RPM); 50000 (4000 RPM) max. Beladung / max. load: 1060 g							
		4214-93	4215	4216	4218	4238			
									
			0519	0521	0523	5127			
						 <sup>12)</sup>			
Kapazität / capacity	ml	4 - 4,5	25	50	100	250			
Maße / dimensions	∅ x L mm	15 x 75	24 x 100	34 x 100	40 x 115	62 x 122			
Anzahl p. Rotor / number p. rotor		120	44	24	16	4			
Drehzahl / speed	RPM	4500	4500	4500	4500	4500			
RZB / RCF	<sup>2)</sup>	4777	4777	4777	4777	4777			
Radius / radius	mm	211	211	211	211	211			
 9 (97%)	sec	115							
 9	sec	116							
Temperatur / temperature	°C <sup>1)</sup>	6							

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















































CP Kulturplatte / Culture plate

DWP Deep Well Platte / Deep well plate



MS Micronic System / Micronic system








4294	4257		4254 + 4255 / 4255-P 8)							
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  90°										
	4259-A 						4449 		4430 	
	Hitachi-Racks 		0554 	0512 	4239 	Corning 	Corning 	Nal-gene 	Nunc 	
			12) 	12) 	12) 					
Kapazität / capacity	ml		650	750	1000	500	250	175	200	
Maße / dimensions Ø x L	mm	20 x 118 x 70	20 x 118 x 70	97 x 139	97 x 152	96 x 176	96 x 147	60 x 162	61,5 x 144,3	60 x 130
Anzahl p. Rotor / number p. rotor		20	20	4	4	4	4	4		4
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500		4500
RZB / RCF	2)	4822	4867	5184	5184	5184	5184	5184		5184
Radius / radius	mm	213	215	229	229	229	229	229		229
 9 (97%)	sec	115								
 9	sec	116								
Temperatur / temperature	°C 1)	7			6					

4294	4254	4254 + 4255 / 4255-P 8)								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  90°										
	4432 		4433 				4434 			
	4255 / 4255-P 8) 		0553 	0578 			0500 	2079 		
	13) 									
Kapazität / capacity	ml	1000	1,5	2,0	5	7	2,7 - 3	4,5 - 5	9	10
Maße / dimensions Ø x L	mm	98 x 138	11 x 38		12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	17 x 70
Anzahl p. Rotor / number p. rotor		4	168		120	120	120	120	76	76
Drehzahl / speed	RPM	4500	4500		4500	4500	4500	4500	4500	4500
RZB / RCF	2)	5184	3600/4686		4618	4618	4618	4618	4777	4777
Radius / radius	mm	229	159/207		204	204	204	204	211	211
 9 (97%)	sec	115								
 9	sec	116								
Temperatur / temperature	°C 1)	6								

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 7) When using removal frame 4259-A please take insert (E2435) out of hanger 4257. Centrifuge removal frame 4259-A only when fully loaded.
- 8) 4255-P: special surface treatment for highest hygienic requirements
- 9) 4255 without lid
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 13) Maximum load 800g. With a load higher than 800g the speed has to be reduced, see label on the bucket. Calculation of the reduced speed see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm<sup>3</sup>".





4294	4254 + 4255 / 4255-P 8)									
<p><b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b></p>  <p>∠ 90°</p>										
	<p>max. Laufzyklen / max. cycles: 30000</p>									
	<p>max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 800 g (4500 RPM), 1000 g (4020 RPM), 1200 g (3670 RPM)</p>									
	<p>4434</p> 									
	0518									
										
Kapazität / capacity	ml	15	10	8	4 - 5,5	9 - 10	4 - 7	8,5 - 10	10	
Maße / dimensions	∅ x L	mm	17 x 100	16 x 80	16 x 81	15 x 75	16 x 92	16 x 75	16 x 100	15 x 102
Anzahl p. Rotor / number p. rotor			76	76	76	76	76	76	76	76
Drehzahl / speed	RPM		4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	<sup>2)</sup>		4777	4777	4777	4777	4777	4777	4777	4777
Radius / radius	mm		211	211	211	211	211	211	211	211
 9 (97%)	sec									115
 9	sec									116
Temperatur / temperature	°C <sup>1)</sup>									6

4294	4254 + 4255 / 4255-P 8)									
<p><b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b></p>  <p>∠ 90°</p>										
	<p>max. Laufzyklen / max. cycles: 30000</p>									
	<p>max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 800 g (4500 RPM), 1000 g (4020 RPM), 1200 g (3670 RPM)</p>									
	<p>4435</p> 		<p>4437</p> 	<p>4438</p> 	<p>4438 + 0726</p> 					
			0509					0519		
										
Kapazität / capacity	ml	2,6 - 3,4	4,9	1,6 - 5	4 - 7	15	25	30	25	
Maße / dimensions	∅ x L	mm	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	25 x 90	25 x 110	24 x 100
Anzahl p. Rotor / number p. rotor			84	84	84	84	48	28	28	28
Drehzahl / speed	RPM		4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	<sup>2)</sup>		4618	4618	4618	4618	4890	4709	4709	4505
Radius / radius	mm		204	204	204	204	216	208	208	199
 9 (97%)	sec									115
 9	sec									116
Temperatur / temperature	°C <sup>1)</sup>									6

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











- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 8) 4255-P: special surface treatment for highest hygienic requirements
- 3) not possible to close the lid

4294		4254 + 4255 / 4255-P 8)						
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  		 						
		max. Laufzyklen / max. cycles: 30000				max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 800 g (4500 RPM), 1000 g (4020 RPM), 1200 g (3670 RPM)		
		4439	4440	4441	4442	4443		
								
		0521	Falcon	Falcon	0513	0526	5127	
			 3)	 3)			 12)	 12)
Kapazität / capacity	ml	50	225	175	50	100	250	290
Maße / dimensions	∅ x L	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 122	62 x 137
Anzahl p. Rotor / number p. rotor		16	4	4	20	8	4	4
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	2)	4573	5184	5184	4890	4551	5003	5003
Radius / radius	mm	202	229	229	216	201	221	221
 9 (97%)	sec						115	
 9	sec						116	
Temperatur / temperature	°C 1)						6	

4274		4275						
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  								
		max. Laufzyklen / max. cycles: 50000 max. Beladung / max. load: 370 g						
		4276-A +0769-A	4277 + 0769-A	0771	4278-A	0703-A		
								
		0531	---	0528				
		 6)	 6)	 6)	 6)			
Kapazität / capacity	ml	100	100	100	50	50		
Maße / dimensions	∅ x L	37 x 200	44 x 168	58 x 161	45 x 130	36,5 x 185		
Anzahl p. Gestell/number p. frame		1	1	1	1	1		
Anzahl p. Rotor / number p. rotor		4	4	4	4	4		
Drehzahl / speed	RPM	2000	2000	2000	2000	2000		
RZB / RCF		961	961	984	984	912		
Radius / radius	mm	215	215	220	220	204		
 9 (97%)	sec						16	
 9	sec						38	
Temperatur / temperature	°C 1)						- 5	

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










CP Kulturplatte /  
Culture plate

DWP Deep Well Platte /  
Deep well plate

MS Micronic System /  
Micronic system





4246													
Winkelrotor 6-fach / Angle rotor 6-times		---		1454		1446		1447		1466		1451	
 <p>45° mit Bioabdichtung / with bio-containment <sup>5)</sup></p>													
		14)	0547	0513	0546	0519	0545	0509					
													
Kapazität / capacity	ml	94	85	50	50	25	30	15	7,5 - 8,5	9 - 10			
Maße / dimensions	Ø x L	mm	38 x 102	38 x 106	29 x 115	29 x 107	24 x 100	26 x 95	17 x 120	15 x 92	16 x 92		
Anzahl p. Rotor / number p. rotor		6		6		6		6		6		6	
Drehzahl / speed	RPM	11500		11500		11500		11500		11500		11500	
RZB / RCF	<sup>2)</sup>	18038		17595		17299		16560		17299		17003	
Radius / radius	mm	122		119		117		112		117		115	
 9 (97%)	sec							64					
 9	sec							64					
Temperatur / temperature	°C <sup>16)</sup>							2					





4246																	
Winkelrotor 6-fach / Angle rotor 6-times		1451		1463		---		1448									
 <p>45° mit Bioabdichtung / with bio-containment <sup>5)</sup></p>																	
		---	---	0518	0521	0548	0549	---									
																	
Kapazität / capacity	ml	8,5 - 10		10		15		50		75		85		10			
Maße / dimensions	Ø x L	mm		16 x 100		15 x 102		17 x 100		34 x 100		35 x 105		38 x 102		16 x 80	
Anzahl p. Rotor / number p. rotor		6		6		6		6		6		6		12			
Drehzahl / speed	RPM	11500		11500		11500		11500		11500		11500		11500			
RZB / RCF	<sup>2)</sup>	17003		17003		17003		17743		17743		18038		17003			
Radius / radius	mm	115		115		115		120		120		122		115			
 9 (97%)	sec							64									
 9	sec							64									
Temperatur / temperature	°C <sup>16)</sup>							2									

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2) Centrifuge containers of glass will not stand RCF values exceeding 4000  
 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".  
 14) Dimensions with cap 38 x 110 mm  
 16) Lowest sample temperature with pre-cooling and maximum speed (only with cooling centrifuges)



4246 Winkelrotor 6-fach / Angle rotor 6-times	1449				4316 11) Winkelrotor 6-fach / Angle rotor 6-times	4317		
	2078	0536	---			0533	0532	
 45° mit Bioabdichtung / with bio-containment 5)						 45°		
Kapazität / capacity ml	1,5	2,0	0,5	3	Kapazität / capacity ml	50	25	
Maße / dimensions Ø x L mm	11 x 38	11 x 38	10,7 x 36	10 x 60	Maße / dimensions Ø x L mm	38 x 148,5	24 x 146,5	
Anzahl p. Rotor / number p. rotor	24		24	24	Anzahl p. Rotor / number p. rotor	6	6	
Drehzahl / speed RPM	11500	11500	11500	11500	Drehzahl / speed RPM	2000	2000	
RZB / RCF 2)	17299	17299	17299	17299	RZB / RCF	805	783	
Radius / radius mm	117	117	117	117	Radius / radius mm	180	175	
 9 (97%) sec	64				 9 (97%) sec	16		
 9 sec	64				 9 sec	29		
Temperatur / temperature °C 16)	2				Temperatur / temperature °C 1)	- 10		

4248 Winkelrotor 14-fach / Angle rotor 14-times	1473				1474	1472	---		
	0518				0509	0519	0545	0513	
 36° mit Bioabdichtung / with bio-containment 5)									
Kapazität / capacity ml	15	10	9 - 10	7,5 x 8,5	8,5 - 10	15	25	30	50
Maße / dimensions Ø x L mm	17 x 100	15 x 102	16 x 92	15 x 92	16 x 100	17 x 120	24 x 100	26 x 95	29 x 115
Anzahl p. Rotor / number p. rotor	14	14	14	14	14	14	14	14	14
Drehzahl / speed RPM	10000	10000	10000	10000	10000	10000	10000	10000	10000
RZB / RCF 2)	16658			16658	16658	16994	16211	17441	17441
Radius / radius mm	149			149	149	152	145	156	156
 9 (97%) sec					103				
 9 sec					106				
Temperatur / temperature °C 16)					6				

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 11) Can only be used in temperatures up to a maximum of 40°C / 104°F.
- 16) Lowest sample temperature with pre-cooling and maximum speed (only with cooling centrifuges)

<b>4266</b> <b>Winkelrotor 6-fach /</b> <b>Angle rotor 6-times</b>  mit Bioabdichtung / with bio-containment <sup>5)</sup>		---		5641	5642	5643	5644		5646		
		5127	---	0545	0519	0546	14)	0547	0549	0518	
Kapazität / capacity	ml	250	10	30	25	50	94	85	85	15	
Maße / dimensions	∅ x L	61 x 122	16 x 80	26 x 95	24 x 100	29 x 107	38 x 102	38 x 106	38 x 102	17 x 100	
Anzahl p. Rotor / number p. rotor		6	48	18		6	6		6	42	
Drehzahl / speed	RPM	9500	9500	9500		9500	9500		9500	9500	
RZB / RCF	<sup>2)</sup>	14025	13420	12915		12108	12310		12310	13319	
Radius / radius	mm	139	133	128		120	122		122	132	
 9 (97%)	sec	82									
 9	sec	96									
Temperatur / temperature	°C <sup>16)</sup>	2									

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2) Centrifuge containers of glass will not stand RCF values exceeding 4000

5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

15) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

14) Dimensions with cap 38 x 110 mm

16) Lowest sample temperature with pre-cooling and maximum speed (only with cooling centrifuges)